

Consistent Quality.



DiaScreen® Reagent Strips for Urinalysis The Quick and Easy Screening Tool for Accurate UAs


- High quality test strips
- Easy to read – test strip spacing matches color chart
- Only two read times for visual use
- Manufactured in the USA
- Large selection of test strip configurations
- Correlates well with leading brands of reagent strips
- May be used visually or with the DiaScreen®50 Urine Chemistry Analyzer



DIASCREEN®

DiaScreen® Liquid Controls


- Convenient and ready to use
- No measuring or mixing required
- Positive and negative controls help comply with quality assurance testing standards
- Compatible with major test strip brands



Cross Reference Guide

Hypoguard	Product #	Bayer Diagnostics	Product #	Roche Diagnostics	Product #
DiaScreen® 10 (glucose, protein, ketone, pH, blood, bilirubin, urobilinogen, nitrite, specific gravity, and leukocytes)	D11000	Multistix® 10 SG	2161	Chemstrip® 10 with SG	417145
DiaScreen® 9 (glucose, protein, ketone, pH, blood, bilirubin, urobilinogen, nitrite, and leukocytes)	D11900	Multistix® 9	2162	Chemstrip® 9	417109
DiaScreen® 8 (glucose, protein, ketone, pH, blood, bilirubin, urobilinogen, nitrite)	D11800	N-Multistix®	2178		
DiaScreen® 7 (glucose, protein, ketone, pH, blood, bilirubin, urobilinogen)	D11700	Multistix®	2179	Chemstrip® 7	417126
DiaScreen® 6 (glucose, protein, ketone, pH, blood, bilirubin)	D11600	Bili-Labstix®	2180	Chemstrip® 6	417135
DiaScreen® 5 (glucose, protein, ketone, pH, blood)	D11500	Labstix®	2181		
DiaScreen® 4ObL (glucose, protein, blood, leukocytes)	D12400			Chemstrip® 4OB	417144
DiaScreen® 4NL (glucose, protein, nitrite, leukocytes)	D13400	Uristix® 4	2166		
DiaScreen® 4pH (glucose, protein, pH, blood)	D11400	Hema-Combistix®	2182		
DiaScreen® 3 (glucose, protein, pH)	D11300	Combistix®	2183		
DiaScreen® 2GK (glucose, ketone)	D12200	Keto-Diastix®	2882	Chemstrip® uGK	00513
DiaScreen® 2GP (glucose, protein)	D11200	Uristix®	2184	Chemstrip® 2GP	200743
DiaScreen® 1K (ketone)	D12100	Ketostix®	2881	Chemstrip® K	00515
DiaScreen® 1G (glucose)	D11100	Diastix®	2803	Chemstrip® uG	00511
DiaScreen® 1B (blood)	D13150	Hemastix®	2190		
DiaScreen® Liquid Controls (2 x 12 mL positive and negative controls)	D76004				
DiaScreen® 50 Urine Chemistry Analyzer	D10000				

Correlation Analysis



The chart above depicts the comparison of DiaScreen® 10 Reagent Strips for Urinalysis to Bayer Multistix® 10 SG and Roche Chemstrip® 10 with SG. In general, DiaScreen correlates very closely to the leading brands in the market.

Multistix is a trademark of Bayer Corporation and Chemstrip is a trademark of a member of the Roche Group.



Minneapolis, MN 55439 USA
 Phone: 800-818-8877
 Fax: 952-646-3110
www.hypoguard.com

P/N 6006-05
 Rev 09/03